



NC TECHNOLOGIES

Innovative Elemental μ -Analysis

UNIVERSAL ORGANIC MICRO ANALYSIS

ECS 80 SERIES

ECS 8060:

Application and examples


ECS 80 SERIES

ECS 8060



ECS 8060 – N configuration

8060
N, PROTEINS



- ✓ Configuration: N
- ✓ Furnaces: No. 2
- ✓ Auto-Sampler: Pneumatic
- ✓ H₂O and CO₂ chemical traps
- ✓ Detector: TCD

ECS 8060 – Technical Features

Type	N – analysis;
Analysis Time	N – from 4 min up to 8 min;
Analytical Range	From 0,002 up to 20 mg;

Accuracy* <0,2% (certified standard; purity >99.9%)

Precision* <0,1% (certified standard; purity >99.9%)

Sample size: Up to 500 mg (depending on sample size and nature)

Sample type: Liquid and Solid

Pneumatic autosampler: up to 147 positions,

Electronic autosampler: 32, 50 100 positions,

Manual sampler

Dual Furnace System,

Safety Quick-fit Integrated,

H₂O trap, CO₂ trap,

Touch-screen display,

Stand-by mode

* Accuracy and precision are related to samples nature and homogeneity.

ECS 8060 – Technical Features

Physical Specifications

Dimensions 110 x 56 x 37 cm,
Weight 90 kg,
Power Supply 230 V,
Absorbed Power 6A,

Gas Requirements:

Helium (99.999% purity), 3-5 bar or
Argon (99.999% purity), 3-5 bar

Oxygen (99.999% purity), 3-5 bar

Air (oil free compressed air)

ECS 8060 – Technical Features

Analytical Conditions

Gas carrier Helium or Argon,

Leak test Automatic leak test,

Furnace temperature Left furnace: 1100°C max

Right furnace: 1100°C max

Oxygen Automatically calculated by the oxygen doser,

Flow Rate Electronic control,

Water trap efficient and easy maintenance,

Detector High Sensitivity TCD,

Software data analysis EAS Clarity,

Calibrations: Linear, Quadratic, cubic

ECS 8060 – Main Application fields



✓ Organic chemistry and pharmaceutical



✓ Soil science and geology- marine science
(distinction between organic and inorganic carbon is available through the previous acidification of the sample)



✓ Environmental analysis



✓ Petrol chemistry and energy



✓ Materials characterization

✓ Food

(Special configuration for big size sample is easily available)


N in food and pet food

The protein content of food is normally determined based on its nitrogen content; this is multiplied by the coefficient 6.25, thus obtaining the crude protein.

In reality, part of the nitrogen of fodder and other foods is not protein in nature, deriving from free amino acids, amides, various organic compounds (nitrogenous bases), as well as from ammonia compounds.


By using ECS8060, nitrogen determination in pet food is easy and low-cost. Thanks to its versatility, both wet and dried samples can be analyzed, giving always high analysis performances and accuracy

ECS 8060 – N, proteins configuration – Typical chromatogram




Standard: Acetalinide

ECS 8060 – N, proteins configuration – Typical calibration curve




Calibration curve for N; range 1mg to 6mg
Standard: Acetanilide


ECS 8060 – Standard example – Sulfanilamide



ECS 8060 – Standard example –Atropine




ECS 8060 – Standard example – Thiourea



ECS 8060 – EDTA vs wheat flour sample

High repeatability and reproducibility of analysis




Summary Table

Sample ID	Sample Weight [mg]	Nitrogen			
		Area [mV.s]	Response	Weight [mg]	Weight [%]
farina_08_01_2018 15_23_22_024	59,672	569,956	569,956	1,023	1,71
Calib\edta 5-6-30 15 ml_08_01_2018 15_14_45_02	56,277	3007,217	3007,217	5,397	9,59

ECS 8060 – EDTA vs dog crunches sample

Possibility of high sample dosage with the same analysis performances




Summary Table

		Sample ID	Sample Weight [mg]	Nitrogen			
				Area [mV.s]	Response	Weight [mg]	Weight [%]
pet food _10_01_2018 16_38_32_024	Colbrick - 1	pet food	158,146	5774,050	5774,050	6,347	4,01
Calib\edta _10_01_2018 15_38_04_019	Colbrick - 1	edta	70,878	6184,005	6184,005	6,797	9,59


ECS 8060 – N, proteins configuration – Meat (wurstel)

	Sample weight	Nitrogen	Nitrogen
Wurstel	[mg]	weight [mg]	weight [%]
N analysis	47.511	1.037	2.18
	71.333	1.85	2.59
	102.883	2.286	2.22
		Average	2.33
		Std. Dev.	0.23
		RDS %	9.70




ECS 8060 – N, proteins configuration – Cheese

	Sample weight	Nitrogen	Nitrogen
Cheese	[mg]	weight [mg]	weight [%]
N analysis	76.58	4.714	6.16
	105.81	7.035	6.65
		Average	6.41
		Std. Dev.	0.35
		RDS %	5.4



ECS 8060 – N, proteins configuration – Dog crunches


Dog crunches	Sample weight [mg]	Nitrogen weight [mg]	Nitrogen weight [%]
N analysis	47.907	2.191	4.57
	75.106	3.703	4.93
		Average	4.75
		Std. Dev.	0.25
		RDS %	5.36



Summary Table						
Sample ID	Sample Weight [mg]	Area [mV.s]	Nitrogen		Element Name	
			Weight [mg]	Weight [%]		
croccantini cane_30_08_2017 11_25_46_023	47,907	6648,827	2,191	4,57	Nitrogen	
croccantini cane_30_08_2017 14_41_17_032	75,106	11238,633	3,703	4,93	Nitrogen	

ECS 8060 – N, proteins configuration – Tomato paste


Tomato Paste	Sample weight [mg]	Nitrogen weight [mg]	Nitrogen weight [%]
N analysis	47.746	0.492	1.03
	80.075	0.796	0.99
		Average	1.01
		Std. Dev.	0.03
		RDS %	2.80



Sample ID	Sample Weight [mg]	Nitrogen			Element Name
		Area [mV.e]	Weight [mg]	Weight [%]	
conc pomodoro_30_08_2017 12_40_54_027 Colibrick - 1	47.746	1492,715	0,492	1,03	Nitrogen
con pomodoro_30_08_2017 14_23_42_031 Colibrick - 1	80,075	2436,473	0,796	0,99	Nitrogen

ECS 8060 – N, proteins configuration – Milk

	Sample weight	Nitrogen	Nitrogen
Milk	[mg]	weight [mg]	weight [%]
N analysis	44.215	0.262	0.59
	77.81	0.487	0.63
		Average	0.61
		Std. Dev.	0.03
		RDS %	4.64





NC TECHNOLOGIES
Innovative Elemental μ -Analysis



Via Milano,15/A - 20041 Bussero (MI), Italy



Phone: +39 02 950 34 69



www.nctechnologies.it